

ONE Therapy Unit Compatible with ALL PREVENA™ Therapy Dressings

DESCRIPTION:

Negative Pressure Options:

- Pre-set, continuous NPWT at -125mmHg for up to 14 days (with dressing changes every 7 days)
- Disposable, single patient use
- Rechargeable battery

SPECIFICATIONS:

Weight with empty canister:

0.64lbs (0.29kg)

Dimensions:

Approx 8.9 x 16.3 x 5.49cm



* Actual size of Therapy Unit

PROTECT Surgical Incisions

Indications for Use:

PREVENA PLUS™ 125 THERAPY UNIT (14 DAY) System is intended to manage the environment of clean closed surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of negative pressure wound therapy.


PREVENA™ Therapy System clinical outcomes:

More than 70+ published clinical studies in a variety of surgical incisions and patient populations have demonstrated that the PREVENA™ Incision Management System is associated with positive clinical outcomes compared to standard post-operative dressings.

PREVENA™ Incision Dressings:

- Contain a skin-friendly wicking interface layer, includes 0.019% ionic silver
- Polyurethane foam bolster specially designed for application to the surgical area
- Therapy should be continuous for a maximum of 7 days per dressing

PREVENA™ System manages surgical incisions by:



Ordering information

Item #	Product Name	Qty
PRE4010	PREVENA PLUS™ 125 Therapy Unit (14 Day)	1
PRE4095	PREVENA PLUS™ 150ml Canister	1
PRE1155US	PREVENA™ PEEL & PLACE™ Dressing – 13cm	5
PRE1055US	PREVENA™ PEEL & PLACE™ Dressing – 20cm	5
PRE3255US	PREVENA PEEL & PLACE™ Dressing – 35cm	5
PRE4055US	PREVENA PLUS™ CUSTOMIZABLE™ Dressing	5

Item #	Product Name	Qty
PRE5055	PREVENA RESTOR ARTHRO•FORM™ Dressing – 33cm x 30cm	5
PRE5155	PREVENA RESTOR ARTHRO•FORM™ Dressing – 46cm x 30cm	5
PRE5255	PREVENA RESTOR™ BELLA•FORM™ Dressing – 21cm x 19cm	5
PRE5355	PREVENA RESTOR™ BELLA•FORM™ Dressing – 24cm x 22cm	5
PRE5455	PREVENA RESTOR™ BELLA•FORM™ Dressing – 29cm x 27cm	5
M6275066/10	V.A.C.® Y-Connector	10

For more information, call 800-275-4524 or visit mykci.com

PREVENA PLUS™ 125 Therapy Unit can be continuously applied for up to 14 days, but the PREVENA™ Dressings have a maximum wear time of 7 days.

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for PREVENA™ Therapy. Please consult the applicable PREVENA™ System Clinician Guide US instructions for use for important limitations, contraindications, warnings and precautions prior to application. Rx only.

Copyright 2019 KCI Licensing, Inc. All rights reserved. Unless otherwise noted, all trademarks designated herein are owned or licensed by KCI Licensing, Inc., Systagenix Wound Management, Ltd., or Crawford Healthcare, Ltd. PRA-PM-US-02125 (12/19)

